

Optimising the Derivatisation of Vitamin D Metabolites Within Clinical Research Using Amplifex™ Diene Reagent

Daniel Blake

AB Sciex, Phoenix House Lakeside Drive, Centre Park, Warrington, Cheshire, WA1 1RX

An increased understanding of the importance of vitamin D has led to a growing interest in the analytical methods used in research for the detection of vitamin D and its metabolites. Here a new derivatisation process using Amplifex™ diene reagent is described, which allows the challenges associated with the measurement of vitamin D and analogues to be overcome.

Introduction

In recent years more researchers have been investigating the role of vitamin D and its metabolites with a view to characterise its metabolic pathways, and to identify and measure its biologically significant metabolites.

The blood levels of 25-hydroxyvitamin D (25-OH-D) are the most appropriate indicator of vitamin D status. The biologically active form of vitamin D, DHVD, is produced when 25-OH-D is hydroxylated and is usually present in the plasma in low pg/ml ranges. Two major forms of vitamin D can be found in the body, these include vitamin D₃, which is endogenously produced, and vitamin D₂, which is only present when it has been given in food or as a supplement [1].

There are several challenges that can make accurate measurement of 25-OH vitamin D difficult. Immunoassays are commonly used in clinical research laboratories for the analysis of 25-OH vitamin D. This is because they are relatively simple to perform and allow a high throughput of samples, however, some immunoassays cannot accurately detect 25-OH vitamin D₂ [2]. Mass spectrometry on the other hand can be used to accurately measure 25-OH vitamin D₂ and 25-OH vitamin D₃, giving separate values for both forms of vitamin D and their metabolites.

Despite the greater sensitivity of mass spectrometry for the detection of 25-OH vitamin D₂ and 25-OH vitamin D₃, there are factors that make accurate measurement of vitamin D difficult. For example, vitamin D in its active form (DHVD₃) is low in abundance (<100 pg/ml) and it has a neutral charge making it difficult to ionise




Figure 1. Fragmentation patterns observed for PTAD or Amplifex™ diene derivatised 1 α , 25-di-OH-vitamin-D (DHVD₃). Amplifex™ diene improved the fragmentation pattern of the analyte.

for mass spectrometry based methods. DHVD₃ is also dehydration prone, and thermally and photochemically labile. The different metabolites of Vitamin D can be similar in structure and molecular weight to DHVD₃ causing interference in the accurate identification and quantification of metabolites by mass spectrometry (or immunoassay). The different isobaric metabolites have to be separated for accurate measurement of DHVD₃ levels. It is also essential to improve the ionisation efficiency and stability of these compounds via derivatisation prior to the mass


spectrometry analysis for a highly selective and reliable detection.

Experimental Methods

The sample preparation, derivatisation and analysis of DHVD₃ can be made simpler, faster and more affordable by using one-step SPE with two inexpensive cartridge and solvents, followed by a one-step derivatisation process.

Traditionally, Cookson-type reagents (PTAD, MBOTAD, DMEQTAD and similar molecules)

2(A)



2(B)




Figure 2. Signal intensity of 1 α , 25-di-OH-vitamin-D (DHVD₃) produced following derivatisation with Amplifex™ diene or PTAD using the following high flow conditions: 2.1 x 100mm, C18, 1.7 μ m, low particle size column e.g. Phenomenex Kinetex or similar, flow rate 0.7 ml/min, column temp 40°C, injection volume 20 μ l, mobile phase A: 0.1 (v/v) formic acid in HPLC water, mobile phase B: 0.1 (v/v) formic acid in acetonitrile. A) DHVD₃ derivatised with Amplifex™ diene compared to underivatized analyte. B) DHVD₃ derivatised with PTAD compared to Amplifex™ diene derivatised analyte.

have been used for the derivatisation of DHVD₃ [3]. These reagents improve the stability and ionisation efficiency of the analyte compared with underivatized samples, however, the fragmentation pattern could be made cleaner and ionisation efficiency improved. A cleaner fragmentation pattern and better ionisation efficiency would increase the accuracy of separation of the different vitamin D analogues, and allow them to be reliably identified and quantified [4].

The new Amplifex™ diene reagent is based on established 1,2,4-Triazoline-based Cookson-type reagents, but the inclusion of a quaternary ammonium group and

other modifications gives an improved mass spectrometry response and more controlled multiple reaction monitoring (MRM) fragmentation, resulting in enhanced analytical detection limits. The Amplifex™ diene reagent can be synthesised in high yields with an ionisation enhancing group, the reagent is also very hydrophilic (which is beneficial for liquid chromatography (LC) methods) but soluble in organic solvent. Amplifex™ diene reacts with vitamin Ds (including vitamin D₃, vitamin D₂ and their respective metabolites) within minutes and unreacted reagents elute very early in reversed phase LC.

A major advantage of the one step

derivatisation process using Amplifex™ diene is that it takes 5 min to perform instead of 30 min, and can be automated and carried out using a 96-well plate. The Amplifex™ diene solution is prepared from dried material, added to the prepared (dried) serum samples (50 μ l per sample) and incubated. After a short incubation the reaction is quenched by the addition of water.

Results

Here we show that the LC-MRM analysis can be performed using a 20 μ l injection, a standard water/acetonitrile (0.1% formic acid) gradient and a 0.7 ml/min flow rate. The QTRAP® 5500 LC/MS/MS system with electrospray ion source in MRM mode was used to quantify vitamin D levels following derivatisation.

Figure 1 displays the fragmentation patterns (measured using the QTRAP 5500 system) of underivatized DHVD₃ compared with PTAD derivatised or Amplifex™ diene derivatised DHVD₃.

A cleaner and more defined fragmentation pattern was observed for DHVD₃ derivatised with the Amplifex™ diene reagent compared with the PTAD reagent, with most of the ion current funnelled to the quantifying transition. The whole analyte was included (structure specific) and was practically free of dehydration. In contrast, DHVD₃-PTAD and its fragment ion underwent dehydration and the quantifying ion carried only a small portion of the analyte. Derivatisation using PTAD would therefore make it difficult to separate isomers that have the Q3 ion in common (for example, DHVD₂). The qualifier ion for PTAD derivatised sample was 298.1 or the water loss at 556.5, showing poor specificity using the PTAD reagent.

Figure 2 (A) demonstrates that the Amplifex™ diene reagent increased the signal intensity for DHVD₃ approximately 180-fold compared with the underivatized analyte. The same concentration of analyte (40 pg) was used for derivatised and underivatized preparations. The transitions for underivatized DHVD₃ were 399.6 – 135.3 compared with transitions of 748.6 – 689.5 after derivatisation with Amplifex™ diene. Figure 2 (B) shows the transitions for DHVD₃-PTAD were closer in range to the underivatized analyte i.e. 574.4 – 314.2 than with Amplifex™ diene. The limits of detection of DHVD₃ and DHVD₂ in stripped serum after derivatisation with Amplifex™ diene have also been determined using the QTRAP 5500 system. The assay requires 200 μ l




Figure 3. Limit of quantification for 1 α , 25-di-OH-vitamin-D (DHVD₃) measured using the QTRAP 5500 following derivatisation with Amplifex™ diene.




Figure 4. Limit of quantification for 1 α , 25-di-OH-vitamin-D (DHVD₂) measured using the QTRAP 5500 following derivatisation with Amplifex™ diene.

of sample, which is quite standard for a clinical laboratory. The limit of quantification for detection of DHVD₃ was approximately 6 pg/ml (Figure 3) whereas the limit of

quantification for DHVD₂ was slightly higher at approximately 10 pg/ml (Figure 4). These detection limits for active vitamin D (DHVD) relate back to the absolute concentration

of the sample and are lower than for more traditional methods using Cookson-type reagents.

Conclusions

The Amplifex™ diene process overcomes many of the difficulties associated with accurate measurement of 1 α , 25-di-OH-vitamin-D₃ such as dehydration, low ionisation efficiency, low structural stability and interference from isobaric isomers. This novel derivatisation reagent can be used to simplify and speed up the chromatography of 1 α , 25-di-OH-vitamin-D₃. Analytes derivatised using Amplifex™ diene have an improved fragmentation pattern, which leads to enhanced sensitivity and lower detection limits compared to more traditional methods. The advantages of using Amplifex™ diene over more traditional derivatisation methods or immunoassays for the detection of 1 α , 25-di-OH-vitamin-D₃ may help to encourage an increase in the use of mass spectrometry methods for this clinical research application.

References

1. Ross, C., Taylor, A., Yaktine, A., Del Valle, H. B. (2011). Institute of Medicine (US) Committee to Review Dietary Reference Intakes for Vitamin D and Calcium. National Academies Press.
2. Glendenning, P., Taranto, M., Noble, J. M., Musk, A. A, Hammond, C., Goldswain, P.R., Fraser, W. D., Vasikarna, S. D. *Annals of Clinical Biochemistry* (2006) 43: 23-30.
3. Shimada, K., Oe, T., Mizuguchi, T. *Analyst* (1991) 116: 1393-1397.
4. Vogeser, M. *The Journal of Steroid Biochemistry and Molecular Biology* (2010) 121: 565-573.

For Research Use Only. Not for use in diagnostic procedures.

RUO-MKT-19-1787-A