

Trace Level Analysis of Aqueous Samples by Dynamic Headspace

by Dan Carrier, Applications Chemist, Anatune Ltd

Dan.carrier@anatune.co.uk

Static Headspace analysis [1] is the most common sample introduction technique for volatile analysis by Gas Chromatography (GC). It works by analysing the vapour above the sample at equilibrium. However, for detection of trace level analytes, this may not be the most suitable method and some form of enrichment may be necessary. For certain analytes, over a 100 fold increase in sensitivity can be achieved with Dynamic Headspace (DHS) analysis when compared to Static Headspace analysis. This article compares DHS with other sample enrichment techniques such as Solid Phase Microextraction (SPME), and Stir Bar Sorptive Extraction (SBSE) using Twister stir bars.

Introduction

In this article, a herbal based liquor and three different orange juices were chosen for enrichment purposes. Both sets of drinks are aqueous based and require some form of enrichment and removal of water for trace level analysis by GC.

Slight differences in the ingredients can drastically change the taste of a drink. The unique smell and taste of different brands of drinks often are due to minor differences in their volatile components. Therefore, it is important to maintain a consistent composition of the drinks to ensure customer satisfaction. Many of the analytes in these drinks are volatile and most are present at trace level. Some of the trace level analytes cannot be detected with Static

Headspace analysis. DHS offers an approach to continuously enrich the analytes making it possible to detect and quantify them by Mass Spectrometry (MS).

This article compares DHS with other sample enrichment techniques such as Solid Phase Microextraction (SPME), and Stir Bar Sorptive Extraction (SBSE) using Twister stir bars.

In DHS, a sample is continuously purged with an inert gas, usually the GC carrier gas, and the volatile compounds are continuously retained onto an adsorptive trap. The trap can then be dried to remove any residual water which may have been collected. Reducing the amount of water is necessary to obtain good chromatography of the analytes which have been enriched. Figure 1 shows a schematic view of the automated

DHS Process.

Firstly, the sample is heated to a required temperature. This allows the sample to form an equilibrium between the liquid and the gas phase in the same manner as Static Headspace analysis. A dual needle then pierces the septa and an inert gas, usually helium, is passed through enabling analytes to be trapped onto the adsorptive trap. Some residual water may have been retained on the trap. Therefore, there is an option to dry purge the trapped sample. If this is required, the dual needle pierces a clean vial and the adsorptive trap is then dried with a set volume of helium at a specified temperature. The trap, containing sample analytes, is transferred into the thermal desorption unit (TDU). A fast temperature ramp is used to desorb the extracted analytes from the TDU onto the Cold Inlet System (CIS) which is set at a cold temperature to focus the analytes of interest. Once the analytes have been focussed into a tight band on the CIS, another fast temperature ramp is used to desorb the extracts onto the GC column.

Full evaporation technique (FET) is used to enhance the detection of volatile analytes by DHS. A low volume aliquot of each drink is placed into an empty headspace vial. The vial is heated to 80°C allowing the analytes in the sample to vaporise whilst leaving the low volatile matrix behind. The FET technique is performed by using a small volume of sample and vaporising the analytes in the headspace vial completely, without having to rely on establishing equilibrium between two




Figure 1: Schematic view of the automated DHS process




Figure 2: Example chromatogram of the herbal based liquor by DHS

phases [2, 3]. The low sample volume used in DHS increases the sensitivity of analytes. This phenomenon may seem counterintuitive, however, this is possible by DHS as the technique effectively enriches most of the analytes.

Twister Stir Bar Sorptive Extraction (SBSE) and Headspace Solid Phase Microextraction (Headspace SPME) are techniques also available to enrich analytes of interest. Twister SBSE is a technique used to enrich analytes onto a magnetic stirrer bar. The most common Twister is coated with a polydimethylsiloxane (PDMS) phase. The Twister bar can be either submerged in an aqueous sample, or can be suspended in the headspace. For the purposes of this article, the Twister is submerged in aqueous samples. The aqueous samples are then placed onto a stirrer plate and stirred for a period of time. Each Twister bar should then be dried with tissue. The enriched analytes are subsequently desorbed in the GC inlet. With Headspace SPME, an SPME fibre is introduced into the headspace and analytes partition into the coating on the fibre. After a fixed time, the fibre is removed and desorbed in the GC inlet. A wide range of fibre coatings are available and this can enable selective extraction of target analytes.

Method

The first series of experiments was designed to determine the differences between the different sampling techniques. A herbal based liquor was chromatographed using Static Headspace analysis, Headspace SPME, Twister SBSE, and DHS using the same GC temperature program and MS conditions. 100 μ l of the herbal based liquor was used

for DHS, SPME, and Twister SBSE; 2 ml of the herbal based liquor was used for Static Headspace; 5-6 replicates were performed to gain understanding of precision using each technique. Estragole and Anethole were also quantified in the herbal based liquor by DHS.

Finally, the three different orange juices were chromatographed using the same DHS method as for the herbal based liquor to identify if differences could be determined between the different sources of juice. 50 μ l of each orange juice was used for DHS; six replicates were performed on orange juice A; individual injections of oranges juices B and C were performed. Table 1 shows the

method conditions used for the analysis.

Results and Discussion

Figure 2 shows an example chromatogram of the herbal based liquor by DHS. Several of the volatile ingredients can be observed by DHS. Only Limonene and Anethole were detected by Static Headspace analysis. Over a 100 fold increase in sensitivity was observed for Anethole using DHS when compared to the Static Headspace analysis for the herbal based liquor. Table 2 shows precision obtained for a number of analytes by DHS, Twister SBSE and Headspace SPME. For DHS, the percentage relative standard deviation (% RSD) for Limonene and Eugenol

Instrumentation	Headspace SPME conditions
Gerstel Multipurpose Sampler MPS 2 XL Gerstel Dynamic Headspace Gerstel Static Headspace Gerstel Thermal Desorption Unit Gerstel Cooled Injection System (CIS) 4 Gerstel SPME Option Maestro Version 1.4.8.14/3.5 Agilent 5975 C inert XL MSD Agilent GC 7890A Twister Stirrer plate	60 °C incubation temperature Fibre: Carboxen/PDMS (incubation time 10 minutes) Inlet temperature 250°C (Splitless)
Column, Oven Program, and MS conditions	Static Headspace conditions
Column: 30 m HP-Innowax (Agilent) di = 0.25 mm df = 0.25 μ m Pneumatics: He, constant flow = 1 mL/min Oven: 40°C (5 min); 10°C/min; 235°C MSD: Scan, 35 - 350 amu	60 °C incubation temperature Injection volume: 2500 μ l CIS 4: Tenax TA liner CIS 4: Temperature Program 40°C; 10°C/s; 300°C (5 min)
Twister SBSE	Dynamic Headspace conditions
100 μ l of the herbal based liquor in 4.9 ml deionized water PDMS Twister (2 cm x 1 mm id) Stirred for 2 hours TDU temperature program 50°C (2 minutes); 720°C/min; 280°C (2 min) CIS 4: Baffled glass liner, CIS 4: Temperature Program -150°C; 16°C/s to 150°C; 12 °C/s to 220°C (5 min)	Trap: Tenax TA DHS Trap temperature 40°C Incubation temperature 80°C Purge volume 50 ml Purge flow 100 ml/min Dry purge volume 950 ml TDU temperature program 50°C; 120°C/min; 350°C (3 min) CIS 4: Tenax TA liner, CIS 4: Temperature Program 10°C; 12°C/s; 250°C (5 min)

Table 1: Instrumentation and method conditions

Technique		Eucalyptol	Isopentylalcohol	Fenchone	Furfural	Linalool	Estragole	Anethole	Eugenol	Limonene
DHS	Mean	123710	1139197	28324	1673222	391793	331255	31165850	6181178	3448051
	sd	6846	67352	1717	166900	30575	22761	885016	1180565	582466
	%rsd	5.5	59	6.1	10.0	7.8	6.9	2.8	19.1	16.9
Twister SBSE	Mean	43033	28265	10382	Not detected	540300	249878	20706991	3812173	Not detected
	sd	3328	1391	1377		37636	15655.2	355162	128999	
	%rsd	7.7	4.9	13.3		7.0	6.3	1.7	3.4	
Headspace SPME	Mean	Not detected	Not detected	Not detected	1911923	220088	437374	25300659	390920	Poor peak shape
	sd				108597	11214	51742	1934005	90742	
	%rsd				5.7	5.1	11.8	7.6	23.2	

%rsd above 10% have been highlighted in red

Table 2: Percentage relative standard deviation (%RSD) for a range of analytes in the herbal based liquor using DHS, Twister SBSE, and Headspace SPME




Figure 3: Linearity for Anethole and Estragole by DHS




Figure 4: Example chromatogram of the herbal based liquor by Headspace SPME showing poor peak shape observed for Limonene

was poor (both >16%). Further method development would be required on these two analytes using DHS to improve the precision. Figure 3 shows the linearity results for Anethole and Estragole. Anethole was present in the herbal based liquor at 13 $\mu\text{g/ml}$ and Estragole was present at 0.16 $\mu\text{g/ml}$. Sensitivity increases were also observed for Twister SBSE and Headspace SPME when compared to Static Headspace analysis. For Headspace SPME, poor peak shape was observed for Limonene. It is thought that this is due to high amounts of water being injected onto the GC. Figure 4 shows the poor peak shape observed for Limonene by Headspace SPME. For Twister SBSE, certain analytes were not detected in comparison to DHS. Furfural was not detected using Twister SBSE. It was thought that Furfural, being hydrophilic, would stay in the highly aqueous solvent and would not interact with the PDMS coating on the Twister bar.

Several of the ingredients in orange juice could be detected by DHS. Figure 5 shows varying levels of the analytes in the three different juices. Six replicas of orange juice A were analysed. Percentage RSDs for Limonene, P-menth-1-en-4-ol, Linalool were all below 5%. Further method development would be required on Eremophilene which was calculated to be 14 % RSD. Precision for β -Pinene was calculated to be 8.5 % RSD.

Conclusion

DHS offers a good automated enrichment process for the volatile components in the herbal based liquor and orange juice. Acceptable quantitation of trace level analytes can be carried out using DHS in aqueous based matrices. Other sample enrichment techniques such as Twister SBSE and Headspace SPME can also be used to give more information on analytes present in

the herbal based liquor and orange juice when compared to Static Headspace analysis.

References

- [1] Static headspace-gas chromatography: theory and practice, Bruno Kolb, Leslie S. Ettre - Wiley-Interscience (2006)
- [2] M. Markelov, J.P. Guzowski, Jr., *Analytica Chimica Acta*, 276 (1993) 235
- [3] M. Markelov, O.A. Bershevits, *Analytica Chimica Acta*, 432 (2001) 213





Figure 5: Comparison different orange juices by DHS

Easy HPLC Method Transitions

Thermo Fisher Scientific, Inc provides invaluable assistance with the HPLC Method Development Calculator to chromatographers looking to optimise methods with the latest column technologies. This interactive range of online tools is hosted on the Thermo Scientific Chromatography Resource Centre and helps calculate initial conditions for gradient and isocratic methods. In addition, an LC column selection tool is provided to assist in the choice of the appropriate column for any application.

With the ability to significantly reduce the time spent in experimentally defining optimal analysis condition, the HPLC Method Development Calculator provides users with quick and easy access to the technical expertise of chromatographers at Thermo Fisher Scientific.

For more information please visit www.thermofisher.com



A new perspective in GC analysis

GCxGC

new!
TOFMS
detection

ZX1:

- Thermal modulation for Multi dimensional GC
- Jet Temperature cools to -189 Deg C using liquid nitrogen
- Modulation range of C1 to C55 for organic components

ZX2:

- Thermal modulation for multi dimensional GC
- Jet temperature cools to -90 Deg C without the need for liquid nitrogen by using Closed Cycle Refrigeration
- Modulation range of C7 to C55 for organic components

Flow modulation of Agilent Technologies:

- Based on micro fluidics technology
- Modulation range C1+

Applications:

- Petrochemistry
- Environment
- Metabolomics
- Nutrients, fragrances and flavours

Cedar Court, Grove Park Business Est. T +44(0) 16 288 220 48
White Waltham, Maidenhead, Berks, info@go-jsb.com
SL6 3LV www.go-jsb.co.uk

Agilent Technologies
Channel Partner

INSPIRATION MEETS INNOVATION!